

National Reference Laboratory for bacteriological contamination of bivalve molluscs

Generic protocol - Enumeration of *Escherichia coli* in
bivalve molluscan shellfish by the colony-count
technique (based on ISO 16649-2)

Note: For official control testing in the UK, either this document or the NRL
generic protocol for enumeration of *Escherichia coli* in bivalve molluscan
shellfish by the most probable number (MPN) technique (based on ISO 16649-
3) must be used

Author(s): James Lowther

Issue number: 1

Date: March 2024



© Crown copyright 2024

This information is licensed under the Open Government Licence v3.0. To view this licence, visit www.nationalarchives.gov.uk/doc/open-government-licence/

This publication is available at www.gov.uk/government/publications

www.cefas.co.uk

Whereas every precaution has been taken in the preparation of this document, Cefas cannot be held responsible for the accuracy of any statement or representation made nor the consequences arising from the use of or alteration to any information contained within. This procedure is intended solely as a general resource for professionals in the field operating within the UK and specialist advice should be obtained if necessary. All references to Cefas must be removed if any alterations are made to this publication.

History of Procedure

Issue	Date	Section	Changes
1	20.03.24		Issue 1

Contents

1. Introduction	6
2. Scope	6
3. Principle	7
4. Safety precautions	7
5. Equipment	7
6. Media	7
7. Microbiological reference cultures for performance testing and controls	8
8. Procedures	8
8.1. Sample transport and receipt	8
8.2. Sample storage	9
8.3. Sample selection	9
8.4. Sample preparation	9
8.4.1. Preparation of BMS using the oyster cracker	9
8.4.2. Preparation of BMS using a shucking knife	10
8.5. Dilution and homogenisation	10
8.5.1. Homogenisation in a blender	10
8.5.2. Homogenisation using a stomacher	11
8.5.3. Preparation of 1 in 20 dilution	11
8.6. Inoculation and incubation	11
8.7. Counting the colony forming units	12
8.8. Test results	12
8.8.1. Maximum countable number of colonies	12
8.8.2. Determination of the reliability of colony counts on parallel plates	12
8.8.3. Determination of the reliability of colony counts at successive dilutions	13
8.8.4. Calculation and reporting of results	14

8.8.5. Use of the Shellfish colony-count calculator	16
9. Uncertainty of test results	17
10. Quality control.....	17
10.1. Proficiency testing.....	17
10.2. Trend monitoring.....	18
10.3. Internal Quality Assurance (IQA)	18
11. References	19
12. Appendices.....	21
12.1. Appendix 1: Sample sizes of shellfish required for <i>E. coli</i> analysis	21
12.2. Appendix 2: Limits of acceptable variability for colony counts of two parallel Petri dishes	23
12.3. Appendix 3: Limits of acceptable ratios for colony counts at successive dilutions	24

1. Introduction

The risks of infectious human diseases acquired from the consumption of bivalve molluscan shellfish (BMS) are internationally recognised. These health hazards are largely due to the phenomenon of filter-feeding whereby BMS concentrate and retain various bacterial and viral pathogens, often derived from sewage contamination of their growing waters. The risks of exposure to infectious agents are compounded by the traditional consumption of raw, or only lightly cooked, BMS. Historically, enteric bacteria, such as faecal coliforms, have been adopted as surrogate indicator organisms to assess the quality of bivalve shellfish flesh, and, consequently, to predict the risk of exposure to enteric pathogenic viruses.

In the UK, the criteria for laying down the microbiological standards for BMS are set out in Retained Regulation (EC) 2017/625 (Anon, 2017a) and Retained Regulation (EC) 2073/2005 (Anon, 2005a) for classification and end product testing respectively. In these regulations *Escherichia coli* is used as an indicator of faecal contamination of shellfish and the reference method for enumeration of *E. coli* is stated as the most probable number (MPN) technique based on ISO 16649-3 (Anon, 2015). These regulations further stipulate that alternative methods may be used if they are validated according to the criteria in ISO 16140.

In 2010, the RIKILT Institute of Food Safety, Wageningen University, the Netherlands undertook a validation study (Jacobs-Reitsma and Pol-Hofstad, 2010) to demonstrate the suitability of a colony-count (pour plate) protocol for enumeration of *E. coli* that it had developed based upon ISO 16649-2 (Anon, 2001) as an alternative method to the MPN. This validation was carried out according to ISO 16140:2003 (Anon, 2003). The validation was certified at the time by MicroVal as an independent certification body and further approved by Cefas as the EURL. A further validation study on the colony-count method using the original study data was carried out by the National Institute for Public Health and the Environment, the Netherlands, in 2021 (Pol-Hofstad and Jacobs-Reitsma, 2021) to reflect updated validation criteria in ISO 16140-2:2016 (Anon, 2016).

The colony-count method is currently listed as one of two accepted alternative methods in the EU Community Guide to the Principles of Good Practice for the Microbiological Classification and Monitoring of Bivalve Mollusc Production and Relaying Areas with regard to Implementing Regulation 2019/627 (Anon, 2021).

2. Scope

This protocol has been produced with reference to ISO 16649-2 (Anon, 2001), relevant parts of other ISO standards referenced in Section 11, and the two validation studies carried out in the Netherlands (Jacobs-Reitsma and Pol-Hofstad, 2010; Pol-Hofstad and Jacobs-Reitsma, 2021). In the context of this test detectable *E. coli* expresses β -glucuronidase activity at $44\pm1^\circ\text{C}$.

Note: This method has only been validated for BMS and is not recommended for analysis of other marine shellfish species e.g. gastropods and echinoderms.

3. Principle

The method used to enumerate *E. coli* in BMS is a pour plate, colony-count technique. Diluted bivalve shellfish homogenate is mixed with TBGA/TBX medium in sterile Petri dishes and incubated at $37\pm1^{\circ}\text{C}$ for 4 ± 0.5 hours (resuscitation step). The inoculated TBGA/TBX plates are then transferred to $44\pm1^{\circ}\text{C}$ for a further 18-24 hours. *E. coli* is confirmed by the presence of blue/green colonies within the medium indicative of β -glucuronidase activity.

Note: Strains of *E. coli* that do not grow at $44\pm1^{\circ}\text{C}$ and those that are β -glucuronidase negative, such as *E. coli* O157 and some other strains of pathogenic *E. coli*, will not be detected by this method.

4. Safety precautions

Standard microbiology safety precautions should be applied throughout. Laboratories should perform a full risk assessment before performing this procedure. Homogenisation of bivalve shellfish should be performed in a Class II safety cabinet to reduce the risk of infection from aerosol inhalation. *E. coli* should be handled in accordance with ACDP category 2 guidelines.

5. Equipment

- Waring blender and jars or Stomacher and stomacher bags
- Class II safety cabinet
- Refrigerator at $3\pm2^{\circ}\text{C}$
- Sterile glassware
- Protective gloves – single use
- Safety gloves – for example a chain mail glove
- Incubator at $37\pm1^{\circ}\text{C}$ and $44\pm1^{\circ}\text{C}$
- Electric top pan balance
- Pipette filler and graduated pipettes or automatic pipettor and pipette tips of a range of sizes e.g. 1ml and 10ml. For handling shellfish homogenates (1 in 2 and 1 in 20 dilutions) open-ended pipettes may be necessary due to high concentrations of particulate matter
- Shucking knife, oyster cracker or other suitable equipment for opening and dissecting shellfish
- Food grade plastic bags
- Steamer, microwave or equivalent apparatus for melting prepared media
- Waterbath at $44\text{--}47^{\circ}\text{C}$
- Sterile 90mm Petri dishes

6. Media

- Peptone salt solution (PSS) (referred to as Maximum Recovery Diluent (MRD) in the UK)
- Tryptone bile glucuronide agar (TBGA/TBX)

7. Microbiological reference cultures for performance testing and controls

The NRL recommends the use of the positive and negative controls throughout the procedure. The strains and criteria included in Table 1 are recommended for use as a minimum.

Table 1: Microbiological reference strains used for control purposes

Media type	Control strain	WDCM ^a	Criteria	Characteristic reaction	Strain choice ^c
Tryptone bile glucuronide agar (TBGA/TBX)	<i>Escherichia coli</i>	00202	Blue to blue-green colonies	β-glucuronidase positive (weak)	2
	<i>Escherichia coli</i>	00012 or 00013	Blue to blue-green colonies	β-glucuronidase positive	1
	<i>Enterococcus faecalis</i>	00009 or 00087	Total inhibition	no growth	1
	<i>Pseudomonas aeruginosa</i>	00025	White to beige colonies	β-glucuronidase negative	1
	<i>Citrobacter freundii</i> ^b	00006			

^a Follow the link (<http://www.phe-culturecollections.org.uk/products/bacteria/WDCMStrains.aspx>) to obtain NCTC number for WDCM reference strains.

^b Following the sub-culturing of *C. freundii* onto TBGA/TBX plates, growth is not always present. We therefore recommend the use of *P. aeruginosa* over *C. freundii*.

^c Reference strains given in ISO 16649-3 for use in performance testing. Strain selection: 1 - Laboratory must select 1 strain from the list provided; 2 - Strain to be used as a minimum.

8. Procedures

8.1. Sample transport and receipt

Samples must be placed in an intact food grade plastic bag and properly packed in a cool box with ice packs – packed in this manner they should reach a temperature between 0°C and 10°C within 4 hours and then maintain this for at least 24 hours¹. For samples where less than 4 hours have elapsed between collection from the production area and receipt at the laboratory, the internal cool box air temperature (or between-bivalve shellfish sample temperature) should be less than the temperature recorded at the time of sampling. Samples from harvesting areas should have been rinsed (but not immersed) and drained at time of sampling. A sample should be regarded as unsatisfactory if on receipt at the laboratory the sample is frozen, the container is leaking, the shellfish are covered in mud or immersed in water or mud/sand.

The sample transport criteria given here are extracted from ISO 6887-3 (Anon, 2020). The use of alternate sample transport criteria may be acceptable, where verification studies have been undertaken and the results of those studies demonstrate that there is no significant effect on the quality of the test results. For Official Controls, it is recommended that verification studies

¹ A temperature data logger may be used to monitor the sample temperature during transit. The data logger must be held in a central position within the cool box and not allowed to come into contact with the ice packs. If the logger has not been activated or is found to be in contact with the ice packs on arrival at the laboratory, then the temperature of the shellfish sample itself should be measured using an appropriately calibrated temperature probe.

supporting the use of sample transport criteria outside of the ranges given in ISO 6887-3 are approved by the relevant Competent Authority (CA).

Note: The UK NRL has carried out studies on the effect of extended storage time and elevated temperatures on certain species of BMS. The report of this work can be accessed using the following web link <https://www.cefasc.co.uk/nrl/information-centre/nrl-reports/>

8.2. Sample storage

Upon receipt in the laboratory the temperature of the samples should be recorded. Samples should be examined immediately - if storage in the laboratory is necessary then samples should be stored at $3\pm 2^{\circ}\text{C}$ and should be processed within 24 hours of collection. If microbiological analysis cannot be initiated within 24 hours of sample collection, an upper limit of 48 hours should be used.

8.3. Sample selection

Choose shellfish that are alive according to the following criteria:

- Reaction or movement of exposed flesh after touching using a sterile shucking knife.
- Shellfish open and close of their own accord.
- A tap on the shell causes closing or movement.
- Tightly closed shellfish.

Discard all dead shellfish and those with obvious signs of damage. Select the appropriate number depending on the species (Appendix 1). More shellfish can be used, if necessary, to produce the required volumes for each analysis.

8.4. Sample preparation

Any mud and sediment adhering to the shellfish should be removed prior to opening the shellfish by rinsing/scrubbing under cold, running tap water of potable quality. Shellfish should not be re-immersed in water as this may cause them to open. Open all selected shellfish as described below ².

8.4.1. Preparation of BMS using the oyster cracker

Sterilise the blade of the oyster cracker before use (see Figure 1 and Figure 2). Place a single animal in a weighing dish on the platform underneath the blade. Lower the lever so that the blade engages the hinge of the shellfish, then fully pull the lever down so that the blade separates the shells. Using a sterilised shucking knife cut the muscle and scrape the meat of both shells into the sterilised container. Transfer any liquor collected in the weighing boat into the sterilised container.

² Alternative suitable equipment can be used to open bivalve shellfish.




Figure 1: Oyster cracker apparatus



Figure 2: Opening an oyster using an oyster cracker

8.4.2. Preparation of BMS using a shucking knife

Open all selected shellfish as described below using a sterilised shucking knife or equivalent and empty meat and liquor into a sterile container. If sterilised by heating allow the knife to cool before using. When opening shellfish, ensure that the hand holding the shellfish is protected with a heavy-duty safety glove to prevent cuts.

8.4.2.1. Oysters and clams

Insert the knife between the two shells towards the hinge end of the animal. Push the knife further into the animal and prise open the upper shell, allowing any liquor to drain into the sterilised container. Push the blade through the animal and sever the muscle attachments by sliding across the animal. Remove the upper shell and scrape the contents of the lower shell into the sterilised container.

8.4.2.2. Mussels and cockles

Insert the knife in between the shells of the animal and separate the shells with a twisting motion of the knife. Collect the liquor from the animal in the sterilised container then cut the muscle between the shells and scrape the contents into the sterilised container.

8.5. Dilution and homogenisation

Weigh the sterilised container and calculate the weight of the contents by subtracting the weight of the pre-weighed container to the nearest gram.

Note: Complete either sections 8.5.1 or 8.5.2.

8.5.1. Homogenisation in a blender

Measure 1ml of sterile PSS per 1g of shellfish (± 2 ml) using a measuring cylinder. Place contents of sterilised container into a 1 litre blender jar³ with the PSS from the measuring cylinder and homogenise at high speed for approximately 1 minute (4 bursts of 15 seconds with at least 5

³ If shellfish are particularly small, it may be necessary to use a smaller blender to achieve a consistent homogenate.

seconds between bursts) in a class two microbiological laminar flow cabinet. Decant the contents (1 in 2 dilution) back into the sterilised container.

8.5.2. Homogenisation using a stomacher

Place the contents of the sterilised container into at least three stomacher bags (e.g. place three stomacher bags inside each other to avoid small pieces of shell from puncturing the bags). Remove excess air from the bags and operate the stomacher for 2-3 minutes. Transfer 50g of homogenised shellfish into another set of 3 stomacher bags, add 50ml of PSS and homogenise for a further 2-3 minutes. After homogenisation decant the contents (1 in 2 dilution) into a new sterilised container.

Note: Results from proficiency testing distributions of cockles have tended to give lower *E. coli* concentrations following stomaching than following blending. NRL advice is therefore to use blending for samples of cockles and other small species such as clams.

8.5.3. Preparation of 1 in 20 dilution

Make a further decimal dilution (1 in 20) using the 1 in 2 dilution as prepared in 8.5.1 or 8.5.2 by adding 1ml of the 1 in 2 dilution to 9ml of PSS.

Note: it may be necessary to use an open-ended pipette to transfer shellfish homogenate due to high concentrations of particulate matter.

8.6. Inoculation and incubation

Melt TBGA/TBX medium using a steamer, microwave or equivalent, then leave to cool in a waterbath at 44-47°C.

Inoculate 5 sterile Petri dishes (90mm diameter) with 2ml of the 1 in 2 dilution. Inoculate 2 sterile Petri dishes (90mm diameter) with 2ml of the 1 in 20 dilution. Use a new sterile pipette for each dilution.

Pour approximately 15-18ml melted and cooled TBGA/TBX medium into each Petri dish. Carefully mix the inoculum with the medium by gentle swirling and allow the mixture to solidify, with the Petri dishes standing on a cool horizontal surface.

Incubate the TBGA/TBX plates inverted for 4 hours at 37±1 °C, followed by 18–24 hours at 44±1 °C.

Note: it may be necessary to use an open-ended pipette to transfer shellfish homogenate due to high concentrations of particulate matter.

Note: The time which elapses between the distribution of the inoculum in a dish and pouring of the medium shall not exceed 15 min.

8.7. Counting the colony forming units

After incubation examine all inoculated TBGA/TBX plates for the presence of typical i.e. blue to blue-green colonies. Record the number of typical colonies on each of the five plates inoculated with the 1 in 2 dilution and each of the two plates inoculated with the 1 in 20 dilution. Atypical colonies should not be included in the count.

Note: in this protocol the maximum countable number of typical colonies for a set of plates at a single dilution is set at a level that corresponds to an average of 150 typical colonies per plate i.e. 750 colonies in total for a set of 5 plates or 300 colonies in total for a set of 2 plates (8.8.1). It is not possible for a set of plates to average ≤ 150 typical colonies and to meet the requirements on the reliability of colony counts of parallel plates (8.8.2) where a single plate within the set contains ≥ 192 typical colonies.

Therefore, where a count of ≥ 192 typical colonies is reached for any single plate, the entire set at the affected dilution can be considered to be above the maximum countable number, and the remaining plates in the set do not need to be counted.

8.8. Test results

Test results are generated from the colony counts for all plates (8.7). The counts may be checked against the maximum countable number and for acceptable reliability, and results calculated and reported as described in sections 8.8.1 to 8.8.4. Alternatively, the Shellfish colony-count calculator can be used (8.8.5). This calculator simultaneously takes into account the maximum countable number, carries out checks for acceptable reliability, calculates the result and shows how the result should be expressed when reporting.

8.8.1. Maximum countable number of colonies

In this protocol the maximum countable number of typical colonies for a set of plates at a single dilution is set at a level that corresponds to an average of 150 typical colonies per plate i.e. 750 colonies in total for a set of 5 plates or 300 colonies in total for a set of 2 plates.

Where the maximum countable number of colonies is exceeded for the set of plates using the 1 in 2 dilution, calculation of results should be carried out using the colony count at the 1 in 20 dilution. However, if the results at the 1 in 20 dilution are rejected due to unacceptable variation between the counts on the two plates (8.8.2), then results should be reported as **>15000 *E. coli* per 100g**.

Where the maximum countable number of colonies is exceeded for the set of plates using the 1 in 20 dilution, then results should be reported as **>150000 *E. coli* per 100g**.

8.8.2. Determination of the reliability of colony counts on parallel plates

To ensure reliable results, the pattern of individual colony counts for all parallel plates within a set at a single dilution should be assessed for acceptable levels of variability using a G-test approach (e.g. Sokal and Rohlf, 2012) as recommended in ISO 7218 (Anon, 2013) and detailed (for 2 parallel plates) in ISO 14461-2 (Anon, 2005b).

For 5 parallel plates (1 in 2 dilution), determine the χ^2 statistic for the set using the following equation:-

$$\chi^2 = \left\{ 2 \left[C_a \cdot \ln \left(\frac{C_a}{(\sum C)/5} \right) + C_b \cdot \ln \left(\frac{C_b}{(\sum C)/5} \right) + C_c \cdot \ln \left(\frac{C_c}{(\sum C)/5} \right) + C_d \cdot \ln \left(\frac{C_d}{(\sum C)/5} \right) + C_e \cdot \ln \left(\frac{C_e}{(\sum C)/5} \right) \right] \right\}$$

where

C_a is the count on the first plate, C_b is the count on the second plate,..... C_e is the count on the fifth plate

$\sum C$ is the sum of the counts on all five plates

The variability between the parallel counts is considered acceptable where χ^2 is less than or equal to 13.277 (the critical value of χ^2 for $p=0.01$ and 4 degrees of freedom).

For 2 parallel plates (1 in 20 dilution), determine the χ^2 statistic for the set using the following equation:-

$$\chi^2 = \left\{ 2 \left[C_a \cdot \ln \left(\frac{C_a}{(\sum C)/2} \right) + C_b \cdot \ln \left(\frac{C_b}{(\sum C)/2} \right) \right] \right\}$$

where

C_a is the count on the first plate and C_b is the count on the second plate

$\sum C$ is the sum of the counts on both plates

The variability between the parallel counts is considered acceptable where χ^2 is less than or equal to 6.635 (the critical value of χ^2 for $p=0.01$ and 1 degree of freedom).

Alternatively, for 2 parallel plates the acceptability of the counts can be determined using the table in Appendix 2.

Where colony counts on a set of parallel plates at one dilution are not acceptable according to the above criteria, then calculation of results should be carried out using the colony counts at the other dilution. Where parallel plate counts at both dilutions are not acceptable then the test result should be considered as **Void**. This should occur in <1% of cases.

Note: if using the colony-count calculator (8.8.5), all steps described in this section are carried out automatically.

8.8.3. Determination of the reliability of colony counts at successive dilutions

To ensure reliable results, the ratio between the total counts (the sum of all plates) at the 1 in 2 and 1 in 20 dilutions should be assessed for acceptability using a G-test approach (e.g. Sokal and Rohlf, 2012) as recommended in ISO 7218 (Anon, 2013) and detailed (for successive dilutions where the same numbers of plates are used) in ISO 14461-2 (Anon, 2005b).

Determine the χ^2 statistic for the ratio of the total counts at the two dilutions using the following equation:-

$$\chi^2 = \left\{ 2 \left[\sum_{upper} \cdot \ln \left(\frac{\sum_{upper}}{10.0 \cdot (\sum_{upper} + \sum_{lower}) / 10.4} \right) + \sum_{lower} \cdot \ln \left(\frac{\sum_{lower}}{0.4 \cdot (\sum_{upper} + \sum_{lower}) / 10.4} \right) \right] \right\}$$

where

\sum_{upper} is the total count of the 5 plates at the 1 in 2 dilution

\sum_{lower} is the total count of the 2 plates at the 1 in 20 dilution

The ratio between the total counts at the two dilutions is considered acceptable where χ^2 is less than or equal to 6.635 (the critical value of χ^2 for $p=0.01$ and 1 degree of freedom).

Alternatively, the acceptability of the ratio of the total counts at the two dilutions can be determined using the table in Appendix 3.

Where the ratio between the total colony counts at the two dilutions is not acceptable according to the above criteria, then the test result should be considered as **Void**. This should occur in 1% of cases.

Note: if using the colony-count calculator (8.8.5), all steps described in this section are carried out automatically.

8.8.4. Calculation and reporting of results

Results are calculated using colony counts for all plates retained after the checks described in 8.8.1, 8.8.2 and 8.8.3. Where there are no issues with exceedance of the maximum countable number (8.8.1) or reliability of the counts at either dilution according to 8.8.2 and 8.8.3, then counts from all plates must be used in the calculation even if e.g. no colonies are present at the 1 in 20 dilution.

Where the sum of colony counts, $\sum C$, of all retained plates is ≥ 10 , calculate the result N (in *E. coli* per 100g) as a weighted mean using the following equations depending on which dilutions are retained:-

If plates at 1 in 2 and 1 in 20 dilutions are retained:

$$N = 100 \cdot \left(\frac{\sum C}{5.2} \right)$$

If plates at 1 in 2 dilution only are retained:

$$N = 100 \cdot \left(\frac{\sum C}{5} \right)$$

If plates at 1 in 20 dilution only are retained:

$$N = 100 \cdot \left(\frac{\sum C}{0.2} \right)$$

Round off the results to two significant figures and report as ***N E. coli* per 100g**.

Example:

Where individual plate counts at the 1 in 2 and 1 in 20 dilutions are (10,15,8,11,14) and (0,1) respectively:-

$$N = 100 \cdot \left(\frac{(10 + 15 + 8 + 11 + 14 + 0 + 1)}{5.2} \right)$$

$$N = 100 \cdot \left(\frac{59}{5.2} \right) = 1134.62$$

Results are rounded off to two significant figures and reported as **1100 *E. coli* per 100g**.

Where the sum of colony counts, $\sum C$, of all retained plates is between 4 and 9, calculate the result N (in *E. coli* per 100g) as a weighted mean using the above equations.

Round off the results to two significant figures and report as **N *E. coli* per 100g estimated**.

Example:

Where individual plate counts at the 1 in 2 and 1 in 20 dilutions are (2,1,2,2,1) and (0,0) respectively:-

$$N = 100 \cdot \left(\frac{(2 + 1 + 2 + 2 + 1 + 0 + 0)}{5.2} \right)$$

$$N = 100 \cdot \left(\frac{8}{5.2} \right) = 153.85$$

Results are rounded off to two significant figures and reported as **150 *E. coli* per 100g estimated**.

Where the sum of colony counts of all retained plates is between 1 and 3, the result should be reported as **present, $<N$ *E. coli* per 100g**, where N is calculated using the above equations, substituting 4 for $\sum C$ as follows:-

If plates at both dilutions are retained report the results as **present, <77 *E. coli* per 100g**

If plates at 1 in 20 dilution only are retained report the results as **present, <2000 *E. coli* per 100g**

Where no colonies are present on the retained plates, the result should be reported as **$<N$ *E. coli* per 100g**, where N is calculated using the above equations, substituting 1 for $\sum C$ as follows:-

If plates at both dilutions are retained report the results as **<19 *E. coli* per 100g**

If plates at 1 in 20 dilution only are retained report the results as **<500 *E. coli* per 100g**

Note: if using the colony-count calculator (8.8.5), all steps described in this section are carried out automatically.

8.8.5. Use of the Shellfish colony-count calculator

The shellfish colony-count calculator is available from the NRL website ([NRL Laboratory Protocols - Cefas](#)).

- Download and open a copy of the colony-count calculator from the website.

Note: only cream-coloured cells can be overwritten, all other cells are locked.

- Each row of the calculator will generate results for a single sample (there are enough rows for 100 samples – example details have been added to row 5 but can be overwritten).
- Add the sample name or unique identifier to column B.
- Add the colony counts for the 5 plates at the 1 in 2 dilution to columns D to H.
- Add the colony counts for the 2 plates at the 1 in 20 dilution to columns J and K.

Note: where the maximum countable number for either set was exceeded and the full count not completed (8.7), the set can be “blocked out” on the calculator by adding any number ≥ 151 to all relevant columns.

- Where the variability between counts on parallel plates within a set is unacceptable, “!!!!” will appear in column Y (1 in 2 dilution) and/or column AG (1 in 20 dilution).
- Where the ratio between total counts at the two dilutions is unacceptable, “!!!!” will appear in column AP.
- Raw and rounded results calculated using the 1 in 2 dilution only, the 1 in 20 dilution only and both dilutions are shown in columns AW to BJ.
- Results (to be reported in ***E. coli* per 100g**), taking into account exceedance of maximum countable numbers at one or both dilutions, all checks for acceptability, expression of results at low colony counts (<10) and rounding of results to two significant figures, are shown in column CD. These result cells are conditionally colour-coded according to the classification bracket in which the results sit:-
 - ≤ 230 – shaded green
 - 240 - 700 – shaded yellow
 - 710 - 4600 – shaded peach
 - 4700 - 46000 – shaded pink
 - >46000 – shaded red

Results that are ambiguous in terms of classification brackets (>15000, <2000, <500) are shaded grey and void results are shaded black with yellow text.

Note: where final results are shown in the form “***N* estimated**”, these should be reported as “***N* *E. coli* per 100g estimated**”.

- Lower and upper 95% confidence limits (section 9) are provided in columns CF and CG for information.

9. Uncertainty of test results

Uncertainty inherent in any test method, i.e. instruments, media, analyst performance etc can be assessed by the repeatability and reproducibility of test results. These should be monitored through control tests analysed alongside sample tests, through in-house comparability testing between analysts and through external inter-comparability exercises, which would highlight any uncertainties within the test methods.

The inclusion of the measurement of uncertainty (MU) or confidence limits for the result is not required when reporting but it is recommended that laboratories determine, as part of their quality procedure, the MU for the reported results and provide this information on request.

ISO 7218 (Anon, 2013) includes a formula for determination of theoretical confidence intervals for colony-count methods; the confidence interval becomes narrower in relative terms as the colony count increases. According to this formula, the confidence interval for a result of 230 *E. coli* per 100g obtained using results at both 1 in 2 and 1 in 20 dilutions and derived from a colony count of 12 is 100 – 360. The confidence interval for a result of 4600 *E. coli* per 100g derived from a colony count of 239 is 4000 – 5200. These intervals are included as an output on the calculation spreadsheet (8.8.5).

ISO 19036 (Anon, 2006) gives guidance on the estimation of MU for quantitative determinations in food microbiology.

10. Quality control

Quality control systems offer a laboratory a mechanism to control preventable inaccuracies within a procedure.

10.1. Proficiency testing

In order to comply with the requirements of Retained Regulation (EC) No. 2017/625 (Anon, 2017a), laboratories undertaking microbiological examination of official control samples of shellfish for *E. coli* and/or *Salmonella* spp. are expected to take part in the UKHSA/Cefas Shellfish EQA scheme (for further information contact foodeqa@ukhsa.gov.uk) and NRL proficiency test distributions⁴. These schemes provide an independent assessment of a laboratory's performance against other

⁴ It was agreed at the Laboratory network meeting in 2018 that laboratories analysing OC samples should take part in at least 2 EQA (UKHSA/Cefas) schemes and 1 whole animal distribution per year.

participants and can help improve the performance of the laboratory.

10.2. Trend monitoring

To ensure continuing intra-comparison of test results, trend analysis should be undertaken by regularly reviewing laboratory performance in the Shellfish EQA scheme. Results should be assessed against the participants' median and plotted graphically showing performance over time. This enables recognition of unusual trends in performance compared to other laboratories and allows for appropriate follow-up action. An Excel spreadsheet for this purpose is given at the Shellfish EQA Scheme web page (<https://www.gov.uk/government/publications/shellfish-scheme-trend-analysis>).

In addition, laboratories should keep a check on the number of results which fail acceptability checks for variability of parallel counts at the 1 in 2 and 1 in 20 dilutions or for the ratio of total counts at successive dilutions. For each of these checks, it is expected that 1% of samples will provide unacceptable results. If a higher proportion is observed, this should be investigated.

10.3. Internal Quality Assurance (IQA)

It is recommended that regular (e.g. monthly) monitoring using known levels of target organism are examined to ensure routine *E. coli* procedures continue to be efficient and effective. An example for assessing quantitative methods is the use of Lenticule™ discs.

11. References

- Anon, 2001. ISO 16649-2:2001. Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of β -glucuronidase-positive *Escherichia coli* – Part 2: Colony-count technique at 44 degrees C using 5-bromo-4-chloro-3-indolyl beta-D-glucuronide. International Organization for Standardization, Geneva.
- Anon, 2003. ISO 16140:2003. Microbiology of food and animal feeding stuffs – Protocol for the validation of alternative methods. International Organization for Standardization, Geneva.
- Anon, 2005a. Retained Regulation (EC) No 2073/2005 of the European Parliament and the Council, 15 November 2005 on microbiological criteria for foodstuffs.
- Anon, 2005b. ISO 14461-2:2005. Milk and milk products - Quality control in microbiological laboratories - Part 2: Determination of the reliability of colony counts of parallel plates and subsequent dilution steps.
- Anon, 2006. ISO/TS 19036:2006 Microbiology of food and animal feeding stuffs -- Guidelines for the estimation of measurement uncertainty for quantitative determinations. International Organization for Standardization: Geneva, Switzerland.
- Anon, 2013. ISO 7218:2007 + Amd 1:2013, Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations (as amended). International Organization for Standardization, Geneva.
- Anon, 2015. ISO 16649-3:2015. Microbiology of the food chain - Horizontal method for the enumeration of β -glucuronidase-positive *Escherichia coli* - Part 3: Detection and most probable number technique using 5-bromo-4-chloro-3-indolyl- β -D-glucuronide. International Organization for Standardization, Geneva.
- Anon, 2016. ISO 16140-2:2016. Microbiology of the food chain – Method validation – Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method. International Organization for Standardization, Geneva.
- Anon, 2017a. Retained Regulation (EC) 2017/625 of the European Parliament and of the Council of 15 March 2017 on official controls and other official activities performed to ensure the application of food and feed law, rules on animal health and welfare, plant health and plant protection products.
- Anon, 2017b. ISO 6887-1:2017. Microbiology of food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination – Part 1: General rules for the preparation of the initial suspension and decimal dilutions. International Organization for Standardization, Geneva.
- Anon, 2020. ISO 6887-3:2017 + Amd 1:2020. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination – Part 3: Specific rules for the preparation of fish and fishery products (as amended). International Organization for Standardization, Geneva.
- Anon, 2021. Community Guide to the Principles of Good Practice for the Microbiological Classification and Monitoring of Bivalve Mollusc Production and Relaying Areas with regard to Implementing Regulation 2019/627, Issue 4. Available online at (aesan.gob.es).
- Jacobs-Reitsma W. F and Pol-Hofstad I. E., 2010. Expert lab report on the MicroVal ISO 16140:2003 validation of the TBX pour plate method (ISO 16649-2) for enumeration of *Escherichia coli* in bivalve molluscs, RIKILT report 2007LR07.

Pol-Hofstad I. E. and Jacobs-Reitsma W. F., 2021. Validation of the TBX pour plate method (ISO 16649-2) for the enumeration of *Escherichia coli* in Live Bivalve Molluscs: renewal study for alignment with EN ISO 16140-2:2016. RIVM report 2021-0127.

Sokal, R. R. and Rohlf, F. J., 2012. Biometry: The Principles and Practice of Statistics in Biological Research (Fourth ed.). New York: Freeman. ISBN-13: 978-0-7167-8604-4.

12. Appendices

12.1. Appendix 1: Sample sizes of shellfish required for *E. coli* analysis

The following sample sizes are recommended for inclusion in the homogenisation step (the recommended number for sampling is 10 % greater to allow for morbidity in a proportion of animals on receipt at the laboratory).

Type	Common name ⁵	Scientific name ⁶	Sample size ⁷
Scallops	Mediterranean scallop	<i>Pecten jacobaeus</i>	10 - 12
	King (Great Atlantic) scallop	<i>Pecten maximus</i>	10 - 12
	Queen scallop	<i>Aequipecten (Chlamys) opercularis</i>	15 - 30
	Variegated scallop	<i>Mimachlamys (Chlamys) varia</i>	10 - 18
Oysters	Pacific oyster	<i>Magallana (Crassostrea) gigas</i>	10 - 18
	Portuguese oyster	<i>Magallana (Crassostrea) angulata</i>	10 - 18
	European flat oyster	<i>Ostrea edulis</i>	10 - 18
Mussels	Blue or common mussel	<i>Mytilus edulis</i>	15 - 30
	Mediterranean mussel	<i>Mytilus galloprovincialis</i>	10 - 30
	Northern horse mussel	<i>Modiolus modiolus</i>	10 - 12
	Bearded horse mussel	<i>Modiolus barbatus</i>	15 - 30
Ark Clams	Ark clam	<i>Barbatia barbata</i>	15 - 25
	Noah's ark shell	<i>Arca noae</i>	15 - 30
Clams	Smooth clam	<i>Callista chione</i>	10 - 30
	Striped venus clam	<i>Chamelea gallina</i>	40 - 80
	Rayed artemis	<i>Dosinia exoleta</i>	18 - 35
	Hard clam	<i>Mercenaria mercenaria</i>	12 - 18
	Palourde (Grooved carpet shell)	<i>Ruditapes (Tapes/Venerupis) decussatus</i>	18 - 35
	Manila clam	<i>Ruditapes (Tapes) philippinarum</i>	18 - 35
	Banded carpet shell	<i>Polititapes (Venerupis) rhomboides</i>	20 - 25
	Pullet carpet shell	<i>Venerupis corrugata (senegalensis)</i>	20 - 25
	Warty venus clam	<i>Venus verrucosa</i>	15 - 30
	Atlantic surf clam (Thick trough shell)	<i>Spisula solida</i>	30 - 50
	Cut trough shell	<i>Spisula subtruncata</i>	70 - 90
	Bean clam	<i>Donax spp.</i>	30 - 50
	Wedge shell clam	<i>Donax trunculus</i>	40 - 80
	Common cockle	<i>Cerastoderma edule</i>	30 - 50
	Soft shell clam (Sand gaper)	<i>Mya arenaria</i>	10 - 12

⁵ Species of relevance provided by European NRLs.

⁶ Scientific names are those used by the World Register of Marine Species at the time of writing, however commonly used alternatives for genus or species names are shown in brackets.

⁷ Number of shellfish to be tested for *E. coli*. Retained Regulation (EC) No. 2073/2005 specifies a minimum of 10 animals of any species should be examined. Sample sizes given in this table have been provided by European NRLs. The weight of shellfish flesh and liquor should be at least 50g for the *E. coli* method (for very small species such as *Donax* spp. a minimum amount of 25g is permitted). For species not given in the table, sufficient shellfish should be opened to achieve this minimum weight of flesh and liquor, with the provision that a minimum of ten animals should be used for very large species. In general, the more shellfish that are included in the initial homogenate, the less the final result will be influenced by the inherent animal-to-animal variation in *E. coli* concentration.

Razor Clams and Geoducks	Pacific geoduck	<i>Panopea generosa</i>	10 - 12
	Sword razor	<i>Ensis ensis</i>	10 - 12
	Razor shell	<i>Ensis magnus (arcuatus)</i>	10 - 12
	Atlantic razor clam (Jackknife clam)	<i>Ensis leei</i>	10 - 12
	Pod razor	<i>Ensis siliqua</i>	10 - 12
	Grooved razor shell	<i>Solen marginatus</i>	10 - 12

12.2. Appendix 2: Limits of acceptable variability for colony counts of two parallel Petri dishes

For a set of two parallel Petri dishes at a single dilution, variability between the plates is acceptable provided that for a given upper plate count as shown in the “Upper colony count” column in the table below, the lower plate count in the set of two plates is no lower than the number shown in the “Lower colony count” column e.g. where the upper count in the set of two plates is 100, the other plate must contain at least 67 colonies for variability to be acceptable.

Colony count		Colony count		Colony count		Colony count	
Upper	Lower	Upper	Lower	Upper	Lower	Upper	Lower
1	0	44	24	87	57	130	92
2	0	45	24	88	58	131	93
3	0	46	25	89	58	132	94
4	0	47	26	90	59	133	95
5	1	48	27	91	60	134	96
6	1	49	27	92	61	135	96
7	1	50	28	93	62	136	97
8	1	51	29	94	62	137	98
9	2	52	29	95	63	138	99
10	2	53	30	96	64	139	100
11	3	54	31	97	65	140	101
12	3	55	32	98	66	141	102
13	4	56	32	99	67	142	102
14	4	57	33	100	67	143	103
15	5	58	34	101	68	144	104
16	5	59	35	102	69	145	105
17	6	60	36	103	70	146	106
18	6	61	36	104	71	147	107
19	7	62	37	105	71	148	107
20	7	63	38	106	72	149	108
21	8	64	39	107	73	150	109
22	9	65	39	108	74	151	110
23	9	66	40	109	75	152	111
24	10	67	41	110	76	153	112
25	11	68	42	111	76	154	113
26	11	69	42	112	77	155	113
27	12	70	43	113	78	156	114
28	12	71	44	114	79	157	115
29	13	72	45	115	80	158	116
30	14	73	46	116	81	159	117
31	14	74	46	117	81	160	118
32	15	75	47	118	82	161	119
33	16	76	48	119	83	162	120
34	16	77	49	120	84	163	120
35	17	78	50	121	85	164	121
36	18	79	50	122	86	165	122
37	19	80	51	123	86	166	123
38	19	81	52	124	87	167	124
39	20	82	53	125	88	168	125
40	21	83	54	126	89	169	125
41	21	84	54	127	90	170	126
42	22	85	55	128	91	171	127
43	23	86	56	129	91	172	128

12.3. Appendix 3: Limits of acceptable ratios for colony counts at successive dilutions

For a set of five parallel Petri dishes at the 1 in 2 dilution and two parallel dishes at the 1 in 20 dilution, the ratio between the total counts at the two dilutions is acceptable provided that for a given total count at the 1 in 2 dilution as shown in the below table, the total count at the 1 in 20 dilution falls within the lower and upper limits in the table e.g. where the total count at the 1 in 2 dilution is 100, the total count at the 1 in 20 dilution must be at least 1 and no more than 10 for the ratio to be acceptable.

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
0	0	1
1	0	1
2	0	1
3	0	2
4	0	2
5	0	2
6	0	2
7	0	2
8	0	2
9	0	2
10	0	3
11	0	3
12	0	3
13	0	3
14	0	3
15	0	3
16	0	3
17	0	3
18	0	4
19	0	4
20	0	4
21	0	4
22	0	4
23	0	4
24	0	4
25	0	4
26	0	4
27	0	4
28	0	4
29	0	5
30	0	5
31	0	5
32	0	5
33	0	5
34	0	5
35	0	5
36	0	5
37	0	5
38	0	5
39	0	5
40	0	6
41	0	6
42	0	6
43	0	6
44	0	6

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
45	0	6
46	0	6
47	0	6
48	0	6
49	0	6
50	0	6
51	0	6
52	0	6
53	0	7
54	0	7
55	0	7
56	0	7
57	0	7
58	0	7
59	0	7
60	0	7
61	0	7
62	0	7
63	0	7
64	0	7
65	0	7
66	0	8
67	0	8
68	0	8
69	0	8
70	0	8
71	0	8
72	0	8
73	0	8
74	0	8
75	0	8
76	0	8
77	0	8
78	0	8
79	0	8
80	0	9
81	0	9
82	0	9
83	0	9
84	0	9
85	1	9
86	1	9
87	1	9
88	1	9
89	1	9

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
90	1	9
91	1	9
92	1	9
93	1	9
94	1	9
95	1	10
96	1	10
97	1	10
98	1	10
99	1	10
100	1	10
101	1	10
102	1	10
103	1	10
104	1	10
105	1	10
106	1	10
107	1	10
108	1	10
109	1	10
110	1	11
111	1	11
112	1	11
113	1	11
114	1	11
115	1	11
116	1	11
117	1	11
118	1	11
119	1	11
120	1	11
121	1	11
122	1	11
123	1	11
124	1	11
125	1	11
126	1	12
127	1	12
128	1	12
129	1	12
130	1	12
131	1	12
132	1	12
133	1	12
134	1	12

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
135	1	12
136	1	12
137	1	12
138	1	12
139	1	12
140	1	12
141	1	13
142	1	13
143	1	13
144	1	13
145	1	13
146	1	13
147	1	13
148	1	13
149	1	13
150	1	13
151	1	13
152	1	13
153	1	13
154	1	13
155	1	13
156	2	13
157	2	13
158	2	14
159	2	14
160	2	14
161	2	14
162	2	14
163	2	14
164	2	14
165	2	14
166	2	14
167	2	14
168	2	14
169	2	14
170	2	14
171	2	14
172	2	14
173	2	14
174	2	15
175	2	15
176	2	15
177	2	15
178	2	15
179	2	15

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
180	2	15
181	2	15
182	2	15
183	2	15
184	2	15
185	2	15
186	2	15
187	2	15
188	2	15
189	2	15
190	2	15
191	2	16
192	2	16
193	2	16
194	2	16
195	2	16
196	2	16
197	2	16
198	2	16
199	2	16
200	2	16
201	2	16
202	2	16
203	2	16
204	2	16
205	2	16
206	2	16
207	3	16
208	3	17
209	3	17
210	3	17
211	3	17
212	3	17
213	3	17
214	3	17
215	3	17
216	3	17
217	3	17
218	3	17
219	3	17
220	3	17
221	3	17
222	3	17
223	3	17
224	3	17

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
225	3	18
226	3	18
227	3	18
228	3	18
229	3	18
230	3	18
231	3	18
232	3	18
233	3	18
234	3	18
235	3	18
236	3	18
237	3	18
238	3	18
239	3	18
240	3	18
241	3	18
242	3	18
243	3	19
244	3	19
245	3	19
246	3	19
247	3	19
248	3	19
249	3	19
250	3	19
251	4	19
252	4	19
253	4	19
254	4	19
255	4	19
256	4	19
257	4	19
258	4	19
259	4	19
260	4	20
261	4	20
262	4	20
263	4	20
264	4	20
265	4	20
266	4	20
267	4	20
268	4	20
269	4	20

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
270	4	20
271	4	20
272	4	20
273	4	20
274	4	20
275	4	20
276	4	20
277	4	20
278	4	21
279	4	21
280	4	21
281	4	21
282	4	21
283	4	21
284	4	21
285	4	21
286	4	21
287	4	21
288	4	21
289	4	21
290	4	21
291	4	21
292	4	21
293	5	21
294	5	21
295	5	21
296	5	22
297	5	22
298	5	22
299	5	22
300	5	22
301	5	22
302	5	22
303	5	22
304	5	22
305	5	22
306	5	22
307	5	22
308	5	22
309	5	22
310	5	22
311	5	22
312	5	22
313	5	22
314	5	23

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
315	5	23
316	5	23
317	5	23
318	5	23
319	5	23
320	5	23
321	5	23
322	5	23
323	5	23
324	5	23
325	5	23
326	5	23
327	5	23
328	5	23
329	5	23
330	5	23
331	5	23
332	5	23
333	6	24
334	6	24
335	6	24
336	6	24
337	6	24
338	6	24
339	6	24
340	6	24
341	6	24
342	6	24
343	6	24
344	6	24
345	6	24
346	6	24
347	6	24
348	6	24
349	6	24
350	6	24
351	6	25
352	6	25
353	6	25
354	6	25
355	6	25
356	6	25
357	6	25
358	6	25
359	6	25

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
360	6	25
361	6	25
362	6	25
363	6	25
364	6	25
365	6	25
366	6	25
367	6	25
368	6	25
369	6	25
370	6	26
371	6	26
372	7	26
373	7	26
374	7	26
375	7	26
376	7	26
377	7	26
378	7	26
379	7	26
380	7	26
381	7	26
382	7	26
383	7	26
384	7	26
385	7	26
386	7	26
387	7	26
388	7	27
389	7	27
390	7	27
391	7	27
392	7	27
393	7	27
394	7	27
395	7	27
396	7	27
397	7	27
398	7	27
399	7	27
400	7	27
401	7	27
402	7	27
403	7	27
404	7	27

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
405	7	27
406	7	27
407	7	28
408	7	28
409	7	28
410	8	28
411	8	28
412	8	28
413	8	28
414	8	28
415	8	28
416	8	28
417	8	28
418	8	28
419	8	28
420	8	28
421	8	28
422	8	28
423	8	28
424	8	28
425	8	28
426	8	29
427	8	29
428	8	29
429	8	29
430	8	29
431	8	29
432	8	29
433	8	29
434	8	29
435	8	29
436	8	29
437	8	29
438	8	29
439	8	29
440	8	29
441	8	29
442	8	29
443	8	29
444	8	29
445	8	30
446	9	30
447	9	30
448	9	30
449	9	30

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
450	9	30
451	9	30
452	9	30
453	9	30
454	9	30
455	9	30
456	9	30
457	9	30
458	9	30
459	9	30
460	9	30
461	9	30
462	9	30
463	9	30
464	9	31
465	9	31
466	9	31
467	9	31
468	9	31
469	9	31
470	9	31
471	9	31
472	9	31
473	9	31
474	9	31
475	9	31
476	9	31
477	9	31
478	9	31
479	9	31
480	9	31
481	9	31
482	9	31
483	10	32
484	10	32
485	10	32
486	10	32
487	10	32
488	10	32
489	10	32
490	10	32
491	10	32
492	10	32
493	10	32
494	10	32

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
495	10	32
496	10	32
497	10	32
498	10	32
499	10	32
500	10	32
501	10	32
502	10	33
503	10	33
504	10	33
505	10	33
506	10	33
507	10	33
508	10	33
509	10	33
510	10	33
511	10	33
512	10	33
513	10	33
514	10	33
515	10	33
516	10	33
517	10	33
518	11	33
519	11	33
520	11	33
521	11	33
522	11	34
523	11	34
524	11	34
525	11	34
526	11	34
527	11	34
528	11	34
529	11	34
530	11	34
531	11	34
532	11	34
533	11	34
534	11	34
535	11	34
536	11	34
537	11	34
538	11	34
539	11	34

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
540	11	34
541	11	35
542	11	35
543	11	35
544	11	35
545	11	35
546	11	35
547	11	35
548	11	35
549	11	35
550	11	35
551	11	35
552	11	35
553	12	35
554	12	35
555	12	35
556	12	35
557	12	35
558	12	35
559	12	35
560	12	35
561	12	36
562	12	36
563	12	36
564	12	36
565	12	36
566	12	36
567	12	36
568	12	36
569	12	36
570	12	36
571	12	36
572	12	36
573	12	36
574	12	36
575	12	36
576	12	36
577	12	36
578	12	36
579	12	36
580	12	37
581	12	37
582	12	37
583	12	37
584	12	37

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
585	12	37
586	12	37
587	12	37
588	13	37
589	13	37
590	13	37
591	13	37
592	13	37
593	13	37
594	13	37
595	13	37
596	13	37
597	13	37
598	13	37
599	13	37
600	13	38
601	13	38
602	13	38
603	13	38
604	13	38
605	13	38
606	13	38
607	13	38
608	13	38
609	13	38
610	13	38
611	13	38
612	13	38
613	13	38
614	13	38
615	13	38
616	13	38
617	13	38
618	13	38
619	13	38
620	13	39
621	13	39
622	14	39
623	14	39
624	14	39
625	14	39
626	14	39
627	14	39
628	14	39
629	14	39

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
630	14	39
631	14	39
632	14	39
633	14	39
634	14	39
635	14	39
636	14	39
637	14	39
638	14	39
639	14	39
640	14	40
641	14	40
642	14	40
643	14	40
644	14	40
645	14	40
646	14	40
647	14	40
648	14	40
649	14	40
650	14	40
651	14	40
652	14	40
653	14	40
654	14	40
655	14	40
656	15	40
657	15	40
658	15	40
659	15	41
660	15	41
661	15	41
662	15	41
663	15	41
664	15	41
665	15	41
666	15	41
667	15	41
668	15	41
669	15	41
670	15	41
671	15	41
672	15	41
673	15	41
674	15	41


1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
675	15	41
676	15	41
677	15	41
678	15	41
679	15	42
680	15	42
681	15	42
682	15	42
683	15	42
684	15	42
685	15	42
686	15	42
687	15	42
688	15	42
689	16	42
690	16	42
691	16	42
692	16	42
693	16	42
694	16	42
695	16	42
696	16	42
697	16	42
698	16	42
699	16	43
700	16	43
701	16	43
702	16	43
703	16	43
704	16	43
705	16	43
706	16	43
707	16	43
708	16	43
709	16	43
710	16	43
711	16	43
712	16	43
713	16	43
714	16	43
715	16	43
716	16	43
717	16	43
718	16	43
719	16	44

1 in 2 Total count	1 in 20	
	Lower limit	Upper limit
720	16	44
721	16	44
722	17	44
723	17	44
724	17	44
725	17	44
726	17	44
727	17	44
728	17	44
729	17	44
730	17	44
731	17	44
732	17	44
733	17	44
734	17	44
735	17	44
736	17	44
737	17	44
738	17	44
739	17	45
740	17	45
741	17	45
742	17	45
743	17	45
744	17	45
745	17	45
746	17	45
747	17	45
748	17	45
749	17	45
750	17	45

Tackling global challenges through innovative science solutions

Cefas, the Centre for Environment, Fisheries, and Aquaculture Science, is an Executive Agency of Defra (the UK Government's Department of Environment, Food and Rural Affairs).

Through innovative solutions and world leading applied science we work to ensure a sustainable future for our rivers, seas and the ocean, supporting healthy and productive marine and freshwater ecosystems.



Pakefield Road, Lowestoft, Suffolk, NR33 0HT

The Nothe, Barrack Road, Weymouth, DT4 8UB

www.cefas.co.uk | +44 (0) 1502 562244



© Crown copyright 2024

